# MEASUREMENT STATIONS HC2-AW-(USB)

## SHORT INSTRUCTION MANUAL

### Probe for humidity and temperature

#### Introduction

The ROTRONIC probe type HC2-AW-(USB) measures the humidity and the temperature of hygroscopic products. The digital probe signal is processed by an indicator such as the HygroLab C1 or the HP23-AW-A, or from a pc (HC2-AW-USB).

The desired scale, 0...100 %rh equilibrium humidity or 0...1 Aw water activity, can be selected through the indicator.


#### Measurement arrangement

For each measurement, a sample cup is necessary.

**Exception:** For the measurement of paper, cardboard and other similar products simply put the probe on the sample.


#### Measurement without temperature control

The whole measuring system, probe + sample cup + sample, adapts to the ambient temperature.

#### Measurement with temperature control

The whole measuring system is temperature stabilized by the connected water circulator.

#### Recommendation

The chosen temperature should not deviate by more than  $\pm 5$  °C from the ambient temperature. Reason: The influence of the ambient temperature on the measuring equipment can not completely be avoided by the temperature control if the temperature difference is too high.

#### Measurement of paper sheets, cardboard sheets etc.


The use of a sample cup is not necessary. Place the probe directly onto the sheet.

#### Starting up HC2-AW-USB

The software HW4 professional with AW-Quick function must be installed properly. For detailed information regarding the installation and handling of the HW4 please refer to the software manual (E-M-HC2-AW-USB-Vxx). Connect the HC2-AW-USB via USB interface to the pc. The pc powers the HC2-AW probe, press the power button to switch the probe on. The red LED starts to illuminate as soon the HC2-AW is powered up and has made connection with the indicator.

#### Starting up HC2-AW

Connect the probe cable with the indicator. The indicator powers the HC2-AW probe, press the power button to switch the probe on. The red LED starts to illuminate as soon the HC2-AW is powered up and has made connection with the indicator.


**Power ON/OFF**  
Press this button to switch the device on or off.

**Operation display**  
As soon as the red LED is lit, the device is ready.

#### General guidelines

To make accurate and reproducible measurement, follow the following general guidelines to measure water activity:


- Measure water activity only in rooms with constant temperature; Ideal are rooms with small daily and yearly temperature variations. Avoid heating of the measuring equipment by sun radiation or radiators, draught, body radiation etc.
- Put the samples into the disposable sample cup, and put the lid on it. The probe and the sample must stay in the same room for acclimation (normally room temperature). A error, which occurs very often in the field is putting the sample a night before the measurement will be done into the fridge without acclimation. A other error is to measure a sample which comes directly from production.
- Use a climate chamber to measure water activity on a other temperature, then the ambient temperature. Set the desired temperature and put the probe and sample into it.
- Avoiding warming up the probe and sample through body heat. Do not hold not to long in your hands.
- Check the proper connection between the probe and the indicator. The red LED illuminates when both, the power and the communication are correct.

#### Measure

- **Put the sample into the sample cup.**

To avoid cross-contamination do not fill more then up to the fill line.

- Put the probe immediately on the sample cup after having filled in the sample in order to avoid humidity exchange with the ambient. The measurement time may vary depends on the probe. Humidity and temperature value are ready to read, when the values are stable for at least 1-2 minutes. Devices with trend indicators will show a stable value with two arrows.

#### Sources of errors

- Temperature influence; for instance unsuitable measuring site; body radiation etc.
- Dirty filter or O-ring
- The accuracy of the probe is outside the admissible limits; in this case the probe has to be checked (calibration / adjustment).
- The probes contain chemicals, which cause a drift or destroy the humidity sensor. The sensor insensitive against chemicals, as long the allowed concentration (continuous operation) will not be exceeded.

Pollutant	Formula	Max.Workplace Concentration (ppm)	Allowed Concentration Continuous Operation		
			(mg/m³)	(ppm)	(mg/m³)
Acetic acid	CH <sub>3</sub> COOH	10	25	800	2000
Acetone	CH <sub>3</sub> COCH <sub>3</sub>	1000	2400	3300	8000
Ammonia	NH <sub>3</sub>	25	18	5500	4000
2-Butanone (MEK)	C <sub>2</sub> H <sub>5</sub> COCH <sub>3</sub>	200	590	3300	8000
Chlorine	Cl <sub>2</sub>	0.5	1.5	0.7	2
Ethanol	C <sub>2</sub> H <sub>5</sub> OH	1000	1900	3500	6000
Ethyl acetate	CH <sub>3</sub> COOC <sub>2</sub> H <sub>5</sub>	400	1400	4000	15000
Ethylene glycol	HOCH <sub>2</sub> CH <sub>2</sub> OH	100	260	1200	3000
Ethylene oxide	C <sub>2</sub> H <sub>4</sub> O	3		3	
Formaldehyde	HCHO	1	1.2	2400	3000
Hydrochloric acid	HCl	5	7	300	500
Hydrogen sulfide	H <sub>2</sub> S	10	15	350	500
Isopropanol	(CH <sub>3</sub> ) <sub>2</sub> CHOH	400	980	4800	12000
Methanol	CH <sub>3</sub> OH	200	260	3500	6000
Nitrogen oxides	NO <sub>x</sub>	5	9	5	9
Ozone	O <sub>3</sub>	0.1	0.2	0.5	1
Petrol		300	1200		150000
Sulfur dioxide	SO <sub>2</sub>	5	13	5	13
Toluene / Xylene	C <sub>6</sub> H <sub>5</sub> CH <sub>3</sub>	100	380	1300	5000

#### Periodic calibration of the probe

Both the Pt100 RTD temperature sensor and the corresponding electronics are very stable and do not normally need to be changed or calibrated after factory calibration. The long term stability of the ROTRONIC Hygromer humidity probes is typically better than 1 %rh per year. For maximum accuracy we recommend calibration of the probe about every 6 to 12 months. More frequent calibration can be necessary in applications where the sensor is exposed to pollutants. The calibration can be performed by the user himself on site or in the laboratory / workshop. For routine calibrations the probe should be checked at least at two points, better at three points. The electronics of the HC2-AW-(USB) do not require calibration in the field. The electronics can not be repaired in the field and should be returned to the manufacturer in the case of problems.

#### Calibration instruction with ROTRONIC humidity standard

Generally the following applies:

- Room temperature: 23 °C  $\pm$ 2K
- Instruments, calibrating device and the humidity standard must be at the same temperature. 1K temperature difference between the humidity standard and sensor result in 6 %rh humidity error
- No draughts in the calibration room
- No direct sunshine
- The values of the humidity standards apply at 23 °C
- If the temperature is different from the standard temperature of 23 °C, the corresponding value of the humidity standard has to be taken into account

#### Recommended calibration values

We recommend to calibrate the HC2-AW-(USB) at 10, 35 and 80 %rh.

#### Preparation

- For the calibration an indicator with UART interface is needed (HP23-AW)
- Make sure the probe is connected to the device. (red LED illuminates)
- Place a textile pad in the disposable sample cup

#### Handling

Hold the ampoule of the selected humidity standard at the neck and bring the entire contents into the further, lower part of the ampoule by swinging or lightly tapping against the glass. Break the ampoule at the marked line and empty it completely on the centre of the textile pad while tapping it lightly. Align the probe so that the opening for the lid points downwards. Screw on the lid with calibration solution and switch on the HC2-AW-(USB).

#### Adaption time

The time to reach the humidity equilibrium depends upon the humidity standard.

0.5 to 65 %rh	45 min.	65 %rh to 94.9 %rh	60 min.	>95 %rh	120 min.
---------------	---------	--------------------	---------	---------	----------

- ! Start to calibrate taking first the lower values and go up.

#### Calibration

The HC2-AW-(USB) is ready for calibration after the adaption time. The calibration procedure varies from indicator to indicator, for proper operation refer to the corresponding indicator manual. After every calibration the textile pad and the disposable cup must be renewed and the lid washed out with water and well dried.

For the calibration with further humidity standards, repeat the calibration steps.

- ! The textile pads are only for onetime using, and can not be used twice.

#### Adjustment

Follow the same procedure as by calibrating the HC2-AW-(USB), only the procedure in the indicator is different. Refer to the corresponding manual.

#### Contamination

Contaminated filters may cause measuring errors and prolong the measuring time. Clean the filter periodically or, if necessary, replace it depending upon the degree of contamination. In order not to damage the sensors, remove the filter for cleaning. Clean the filter with soap water, alcohol or a cleaning agent suitable for removing the contamination and rinse it thoroughly at the end with water. Do not install the filter until it is completely dry. Replace filters which can no longer be cleaned by new ones. Should the sensors be severely contaminated despite of the protecting filter, then we recommend having these replaced by our service department.

#### Accessories and spare parts

Accessories: Disposable sample cup made out of polystrol:

- For cup holder with 14mm fill level  
Order code PS-14 (One bag contains 100pcs)
- For cup holder with 40mm fill level  
Order code PS-40 (One bag contains 100pcs)

Spare parts: Spare filter including safety ring: order code ET-W24-Set

#### Technical data

Operating ranges: Humidity 0...1 aw / 0...100 %rh  
Temperature -40...85 °C

Output signal:  
HC2-AW: Digital UART  
HC2-AW-USB: USB  
Accuracy at 23 °C: Humidity  $\pm$ 0.008 aw /  $\pm$ 0.8 %rh /  
Temperature  $\pm$ 0.1K  
Long-term stability: < 1 %rh / year (Normal conditions)  
Factory adjustment: Adjustment at 23 °C, 10, 35 and 80 rh%  
Connector: Standard HC2-connector

Cable length:  
HC2-AW: 1 m  
HC2-AW-USB: 3 m  
Filter: wire mesh  
Weight: Approx. 540 gr.  
Dimensions: 60mm /  $\varnothing$  68mm

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#### ROTRONIC AG, CH-8303 Bassersdorf
#### ROTRONIC Messgeräte GmbH, D-76275 Ettlingen
#### ROTRONIC Instruments (UK) Ltd, Crompton Fields,
#### ROTRONIC Instrument Corp, NY 11788, USA
#### ROTRONIC South East Asia Pte Ltd, Singapore 339156