

Set of Events in Audit Trail

Below is an overview of the most common events that appear within the Rotronic Monitoring System Audit Trail. This explains the meaning of the event, the possible root causes and proposed actions.

If an event occurs within your system that is not listed, or the event can't be explained by the list below, then please export an unfiltered overview of all events before, during and eventually after the event occurred. Send the event overview to rms@rotronic.ch

Event name	Adjust
Meaning	An adjustment was carried out.
Possible root cause	A user did an adjustment on a measuring point.
Action	No reaction required in case this message is expected.

Event name	Alarm script
Meaning	An alarm script has been run.
Possible root cause	The conditions of an alarm script were fulfilled and the alarm script was executed.
Action	React according to internal SOP's.

Event name	Allocate gateway
Meaning	A different gateway has been allocated to a wireless device.
Possible root cause	1) A new gateway with a better signal than the original one has been detected and therefore chosen by the wireless device 2) The original gateway is offline, and communication changes to another gateway
Action	No reaction required in case this message is expected.
Comment	This event can only occur if the gateways are setup to create a redundancy (same network ID and channel)

Event name	Calibrate
Meaning	A calibration was carried out.
Possible root cause	A user did a calibration on a measuring point.
Action	No reaction required in case this message is expected.

Event name	Calibration reminder
Meaning	A calibration reminder was triggered.
Possible root cause	A user set a reminder for a calibration with a defined time lapse. The lapse was achieved.
Action	1) Calibrate the measuring point in question. 2) Afterwards close the reminder eventually set a new reminder.

Event name	Create element
Meaning	An element has been created.
Possible root cause	A user has created a new element (new user, group, device, measuring point or report).
Action	No reaction required in case this message is expected

Event name	Data gap	
Meaning	The RMS server tries to read out (backfill) data from device after recovering from communication interruption. But the backfilling is not successful.	
Case 1	Behavior	Three "read out data" trials with the "timeout" message. Afterwards the "data gap" event occurs
	Possible root cause	A communication problem between the device and the server during "read out data": 1) Faulty power supply 2) A network issue 3) A gateway issue
	Action	Check the device status and make sure it's not in the device timeout status. Retry to read out the data manually: Event page>Data gap event>Try again
Case 2	Behavior	The "Read out data" trial with "no data" message followed by a "data gap" event.
	Possible root cause	The device has no measured data for the asked time period: 1) Faulty power supply (Battery, external supply, PoE) 2) The device was switched off 3) There was a deviation of the time synchronization
	Action	1) Check in the audit trail for related information such as "measuring time not valid" during the period of the data gap 2) Check in the IT system, whether the server time made a big jump

Event name	Delete element
Meaning	An element has been deleted.
Possible root cause	A user has deleted an element (user, group, device, measuring point or template).
Action	No reaction required in case this message is expected

Event name	Devices timeout	
Meaning	Measured data was not received on the server side for longer than defined timeout limit.	
Case 1	Behavior	Timeout for a certain period
	Possible root cause	For wired devices: 1) A problem with the device power supply 2) A network problem between the server and the device For wireless devices: 1) A problem with device power supply 2) A network problem between the server and the gateway 3) A disturbance in the signal or a weak signal to the wireless devices
	Action	1) Check to see if there are a any validation or maintenance work on-going, react according to the internal SOP's and check that the data was backfilled 2) Check the power supply and network of the gateway and/or devices 3) Check the signal strength of the wireless devices
	Comment	If various wireless networks are running, make sure that the networks are not in conflict with each other.
Case 2	Behavior	Timeout remains
	Possible root cause	1) The battery capacity is too low 2) No power supply 3) No network connection to the server 4) A device malfunction
	Action	Check the device status (hardware, power supply, network connection...).

Event name	Edit element
Meaning	An element has been edited
Possible root cause	A user has edited the property of an element (user, device, measuring point).
Action	No reaction required in case this message is expected.

Event name	Inhibit alarm notification
Meaning	Alarm notifications have been inhibited.
Possible root cause	A user has inhibited an alarm or a warning for a certain time.

Action	No reaction required in case this message is expected.
--------	--

Event name	Inhibit technical alarms
Meaning	Technical alarms have been inhibited.
Possible root cause	A user has inhibited technical alarms (sensor error, devices timeout, low battery) for a certain time.
Action	No reaction required in case this message is expected.

Event name	Log in
Meaning	A user has logged in.
Possible root cause	A user has logged in.
Action	1) No reaction required in case this message is expected. 2) In case of an unexpected "Log in", contact your IT responsible.

Event name	Log out
Meaning	A user has logged out.
Possible root cause	A user has logged out.
Action	1) No reaction required in case this message is expected. 2) In case of an unexpected "Log out", contact your IT responsible.

Event name	Low battery
Meaning	A device battery is low.
Possible root cause	1) The battery is used and must be changed. 2) During a battery change, after the old battery was removed, the new battery was inserted quickly. The capacitor on the device was not discharged completely and therefore the counter was not reset.
Action	Change the battery. After removing the old battery, wait minimum 1 minute or press the button on the device before putting in the new one.
Comment	1) Should the time out event occur frequently, the device battery will decrease at an increased rate due to the fact that the device will try to reconnect to the system. 2) The battery monitoring function of RMS is done with a counter and not an actual measurement of the battery capacity.

Event name	Measured value out of limit
Meaning	The measured value was out of defined warning/alarming limits.
Possible root cause	1) Measuring value too high or too low. 2) The wrong alarm/warning levels have been set.
Action	1) Check measured data during alarm/warning.

	<p>2) In case of suspicion, check alarm/warning levels, delays and hysteresis.</p> <p>3) Acknowledge/Inhibit the warning/alarm according to internal SOP's.</p>
--	---

Event name	Measuring time not valid	
Meaning	The time stamp received from the device (device time) is different to the system time (server time).	
Case 1	Behavior	Only a single message. The difference between the device time and the server time < 1 minute
	Possible root cause	A delay on the measuring data transmission. Service restart on server, server reboot, network issues.
	Action	No reaction required, if this happens seldom and the system works further without any issues. In case the message appears too frequently: 1) Check the settings such as the service restart interval 2) Check the IT infrastructure
Case 2	Behavior	Only a single message. The difference between the device time and the server time > 1 minute
	Possible root cause	The device has restarted and it was using the last known time stamp, which was stored internally
	Action	No reaction required, if this happened after re-powering or re-pairing a device. In case message appears too frequently: 1) Check the IIS settings about the application recycle 2) Check the IT infrastructure
Case 3	Behavior	Consecutive messages every interval
	Possible root cause	Device control crashed due to low power supply: 1) Low battery status. 2) Too low supply voltage in case of an external power supply. 3) The device is subjected to very low temperatures. 4) A device malfunction.
	Action	1) Change the device battery. 2) Remove the device from the very low temperature environment and consider the use of an extension cable. 3) Check the external power supply.

Case 4	Behavior	The message appears for different devices at the same time and the device time on different devices are the same.
	Possible root cause	The server time has been updated. The device time and server time will mismatch.
	Action	Check the IT infrastructure. The server time shall not have big jumps

Event name	Probe replaced
Meaning	A probe was replaced.
Possible root cause	A change in probe serial number was detected.
Action	1) No reaction required in case this message is expected. 2) Acknowledge the event according to internal SOP's.
comment	The event is frequently seen during calibration should the probe be removed from the logger for calibration.

Event name	Read out data (not followed by "Data gap")	
Meaning	RMS detects missing data for a period and tries to read back the missing data from the device.	
Case 1	Behavior	With single devices and seldom occurs.
	Possible root cause	Short-term failure of the wired or wireless connection.
	Action	No reaction required in case this message is expected.
Case 2	Behavior	With single devices and frequently occurs
	Possible root cause	1) The connection between the server and the wired device is disturbed, network issue. 2) A disturbance or weak signal for the wireless devices.
	Action	A frequent "read out data" message is an indication of a network issue (frequent communication interruption). Check the network and/or wireless connection.
Case 3	Behavior	Happens to all measuring points at the same time.
	Possible root cause	1) The server was offline due to an application recycle. It starts to backfill the data when it's online again. 2) The connection to the server was disturbed, network issue.

	Action	<p>1) No reaction required if an expected offline of the server was expected.</p> <p>2) Check the IIS settings about the application recycle.</p> <p>3) Check the network connection.</p>
--	--------	---

Event name	Send e-mail	
Meaning	An e-mail was sent to a certain user due to reminder/alarm/warning/error.	
Possible root cause	A reminder/alarm/warning/error occurred and an e-mail was sent according to the notification settings.	
Action	React based on the e-mail content and/or internal SOP's.	

Event name	Send SMS	
Meaning	An SMS was sent to a certain user due to reminder/alarm/warning/error.	
Possible root cause	A reminder/alarm/warning/error occurred and an SMS was sent according to the notification settings.	
Action	React based on the SMS content and/or internal SOP's.	

Event name	Sensor error	
Meaning	An error has been detected on sensor/probe.	
Possible root cause	<p>1) The sensor/probe is defect.</p> <p>2) The sensor/probe is disconnected or has a bad connection to the device.</p> <p>3) A device malfunction has occurred.</p>	
Action	<p>1) Check the audit trail for other related events.</p> <p>2) Check the physical connection of sensor/probe.</p> <p>3) Try with another sensor/probe if possible.</p>	

Event name	Simulator connected	
Meaning	A device is set to simulator mode.	
Possible root cause	<p>1) For the mini modules, the checkbox for simulator can be selected under Setup>device>device properties.</p> <p>2) For the digital probes, the devices are set to simulator mode via the RMS-Config tool.</p>	
Action	No reaction required in case this message is expected.	

Event name	Software exception	
Meaning	A software exception has occurred.	

Possible root cause	1) After an RMS service restart, several "software exception" messages can be triggered. RMS service restart is controlled via the IIS settings 2) Another software problem
Action	If the "software exception" only occurs seldom and right after the RMS service restart, no reactions are required.

Event name	System
Meaning	An event, that is triggered by the system.
Possible root cause	For example an RMS server stop/restart.
Action	No reaction required in case this message is expected

Event name	Telephone call
Meaning	A telephone call was made to a certain user due to reminder/alarm/warning/error.
Possible root cause	A reminder/alarm/warning/error occurred and a telephone call was made according to the notification settings.
Action	React based on the telephone call content and/or internal SOP's.

Event name	Turn switch off
Meaning	A switch output has been turned off.
Possible root cause	1) The conditions of an alarm script were fulfilled and the alarm script was executed and the switch output was turned off. 2) A user turned off the switch output manually using a virtual button.
Action	No reaction required in case this message is expected.

Event name	Turn switch on
Meaning	A switch output has been turned on.
Possible root cause	1) The conditions of an alarm script were fulfilled and the alarm script was executed and the switch output was turned on. 2) A user turned on the switch output manually using a virtual button.
Action	No reaction required in case this message is expected.