

## Product: HP-GTS

### Hard-and Firmware changes

#### Version 1.1 (Valid from January 2021)

This hard- and software update 1.0 extends all functions and features of the last generation of HP-GTS (Version 1.0) with the following modifications and additions:

#### Enhancements/Improvements:

- New display key pad with different On / Off key
- More durable key pad

#### Bugs fixed:

- None

#### Known bugs:

- None

#### Notes:

- Version 1.1 has the same firmware as Version 1.0.  
No Software upgrade possible from Version 1.0 up to Version 1.1,  
as V1.1 include only hardware changes.

