

January 2015, Bassersdorf, Switzerland

Hardware and firmware compatibility of AirChip3000 products

Dear partner

Please find all information regarding AirChip3000 products and its hardware and firmware compatibility in this document.

Products from factory

All AirChip3000 products which are equipped by factory with new hardware (serial number starts with 2xxxxx) have a firmware v3.2 installed. Currently this is the latest firmware version. Therefore no firmware update is provided via the website.

Firmware downgrade

The firmware of these devices can't be downgraded to older versions.

Devices with older hardware

AirChip3000 products with older hardware can't be updated to v3.x. Firmware version v2.0 is the latest for these products.

Confirmed by the manufacturer



Philipp Storrer

Technical Product Manager