

Product: HP31

Hard-and Firmware changes

Version 1.1 (Valid from September 2020)

This hard- and software update 1.1 extends all functions and features of the last generation of HP31 (Version 1.0) with the following modifications and additions:

Enhancements/Improvements:

- New display key pad with different On / Off key
- More durable key pad

Bugs fixed:

- None

Known bugs:

- None

Notes:

- Version 1.1 has the same firmware as Version 1.0.
No Software upgrade possible from Version 1.0 up to Version 1.1, as V1.1 include only hardware changes.

