

## PROJECT FINANCING

*For regulated and non regulated companies, environmental monitoring is crucial in order to comply with regulations and guidelines to ensure complete product traceability and accountability.*



Many different monitoring solutions exist in the market today, the most cost effective (in terms of initial investment) remains simple data loggers such as the Rotronic HL-1D (<100USD). The required values are taken from the display at regular intervals or can be downloaded from the data logger using the HW4 software (FDA CFR 21 Part 11 compliant). Reporting can only be carried out retroactively.

However, such a solution requires a physical individual to carry out the various tasks, adding various sources of error and increasing the overall risk, leading to a much higher total system cost.

In order to reduce risk to a maximum and decrease the total system cost (mainly due to the physical individual), an automated computerised system is preferred. The Rotronic Monitoring System is a fully automated continuous environmental monitoring system. The various data loggers automatically transmit their data to a server (on premise or SaaS). Reporting and alarming is fully automated, ensuring that the various users are always up to date with the latest information and informed in real time about any discrepancies. The downside of such a system is the initial investment, as such Rotronic offers you the following financing solutions:

One off  
payment

Monthly  
payment scheme

Yearly  
payment scheme